

# Specifications: M260



**SHENZHEN COMEN MEDICAL INSTRUMENTS CO.,LTD**

Floor 10, Floor 11 and Section C of Floor 12 of Building 1A &  
Floor 1 to Floor 5 of Building 2, FIYTA Timepiece Building, Nanhuan  
Avenue, Matian Sub-district, Guangming District, Shenzhen,  
Guangdong, 518106,P.R. China

Tel: +86-755-26408879

Fax: +86-755-26431232

Email: [info@szcomen.com](mailto:info@szcomen.com)

Web: [www.comen.com](http://www.comen.com)

# Syringe Pump

## M260



### Physical Characteristics

Dimensions:	<312mm×128mm×133mm (L*W*H)
Weight	≤1.6kg (Including battery)
Screen Size:	3" Segmented LCD Display, 240×128 pixels
Classification:	Type CF, Defibrillation-proof, Class I, EN1789 Certified

### Operating Conditions

Temperature:	5~40°C
Humidity:	15%~95% (non-condensing)
Pressure:	57.0~106.0kPa

### Storage Conditions

Temperature:	-20~+55°C
Humidity:	10%~95% (non-condensing)
Pressure:	57.0~106.0kPa

### Power

Power requirement:	100-240V~, 50/60Hz±1Hz
Applied Input Power:	0.5A~0.3A
Battery Type:	Rechargeable Lithium-ion battery
Battery Capacity:	Standard 7.2V — 2550mAh Optional 7.2V — 5100mAh
Battery Charging Time:	2550mAh: Not more than 4 hours to full charge when powered off 5100mAh: Not more than 7.5 hours to full charge when powered off
Battery performance:	2550mAh: Up to 6 hours of continuous working (at 5ml/h) 5100mAh: Up to 12 hours of continuous working (at 5ml/h)

### Functional

Working modes:	Rate Mode, Time Mode, Weight Mode, Sequential mode, First dose mode, Ramp mode, Micro mode, Intermittent Mode, Dose Time Mode
----------------	---

### Indicators

Hardware Indicators:	Battery status indicator AC Power indicator Alarm Light bar indicator Power Button backlight
On-Screen Indicators:	Infusion related information: Selected drug, Mode, Left Vol, Volume, Rest Time, Infusion Rate, Alarm status, Battery status, Syringe brand and size, Running (Infusion Status) Bar, WIFI Icon Dynamic Pressure System (DPS)

### Interface

AC power interface
Multifunction Port (RS232, DC Power, Nurse Call)

### Dynamic Pressure System (DPS)

Display:	DPS indicator displayed as visual pressure value for occlusion alarm.
Selectable Levels:	4 Selectable Levels: 150, 300, 525, 900 Unit: mmHg

### Ergonomics

Mounting Solutions:	Allows for mounting through a pole clamp for either stationary use on a bed rail or portability through the use of a trolley (TC800).
Portability:	Comes with a Carrying Handle.

**Working Modes:****Standard Modes:**

Rate Mode, Time Mode, Weight Mode

**Optional Modes:**

Micro Mode, Sequential Mode, First Dose Mode, Ramp Mode, Dose Time Mode, Intermittent Mode

**Rate Mode setting ranges:**

Parameters:	Rate, Time and Volume to Be Infused (VTBI) can be set to automatically calculate the remaining parameter.
Rate range:	0.10~2000.00ml/h;
Rate steps:	At 0.10 ml/h~99.9ml/h, step: 0.01 ml/h At 100.0 ml/h~999.9ml/h, step: 0.1 ml/h At 1000 ml/h~2000 ml/h, step: 1 ml/h
VTBI range:	0.10~9999ml
VTBI steps:	At 0.10 ml~99.9ml, step: 0.01 ml At 100.0 ml~999.9ml, step: 0.1 ml At 1000 ml~9999 ml, step: 1 ml
Time range:	00:00:01~99:59:59 (h:m:s); Step: 1s

**Time Mode setting ranges:**

Parameters:	Rate, Time and Volume To Be Infused (VTBI) can be set to automatically calculate the remaining parameter.
Rate range:	0.10~2000.00ml/h;
Rate steps:	At 0.10 ml/h~99.9ml/h, step: 0.01 ml/h At 100.0 ml/h~999.9ml/h, step: 0.1 ml/h At 1000 ml/h~2000 ml/h, step: 1 ml/h
VTBI range:	0.10~9999ml
VTBI steps:	At 0.10 ml~99.9ml, step: 0.01 ml At 100.0 ml~999.9ml, step: 0.1 ml At 1000 ml~9999 ml, step: 1 ml
Time range:	00:00:01~99:59:59 (h:m:s); Step: 1s

**Weight Mode setting ranges:**

Parameters:	Drug Weight, Drug Volume, Patient Weight, Dose Rate and Volume to Be Infused (VTBI) can be set to automatically calculate the infusion rate.
Rate range:	0.10~2000.00ml/h;
Rate steps:	At 0.10 ml/h~99.9ml/h, step: 0.01 ml/h At 100.0 ml/h~999.9ml/h, step: 0.1 ml/h At 1000 ml/h~2000 ml/h, step: 1 ml/h
VTBI range:	0.10~9999ml
VTBI steps:	At 0.10 ml~99.9ml, step: 0.01 ml At 100.0 ml~999.9ml, step: 0.1 ml At 1000 ml~9999 ml, step: 1 ml
Time range:	00:00:01~99:59:59 (h:m:s); Step:1s

**Sequential Mode setting ranges:**

Parameters:	Time, Rate and Volume to be infused (VTBI) can be set for each sequence to automatically calculate the remaining parameter and display a total VTBI.
Sequences:	Up to 10 sequences per pump.
Rate range:	0.10~2000.00ml/h;
Rate steps:	At 0.10 ml/h~99.9ml/h, step: 0.01 ml/h At 100.0 ml/h~999.9ml/h, step: 0.1 ml/h At 1000 ml/h~2000 ml/h, step: 1 ml/h
VTBI range:	0.10~9999ml. The total sum of VTBI for all sequences cannot exceed 9999ml.
VTBI steps:	At 0.10 ml~99.9ml, step: 0.01 ml At 100.0 ml~999.9ml, step: 0.1 ml At 1000 ml~9999 ml, step: 1 ml
Time range:	00:00:01~99:59:59(h: m: s); Step: 1s

**First Dose mode setting range:**

Parameters:	Main time, Main rate, Main volume, First time, First rate and First volume can be set.
Rate range:	0.10~2000.00ml/h;
Rate steps:	At 0.10 ml/h~99.9ml/h, step: 0.01 ml/h At 100.0 ml/h~999.9ml/h, step: 0.1 ml/h At 1000 ml/h~2000 ml/h, step: 1 ml/h
VTBI range:	0.10~9999ml
VTBI steps:	At 0.10 ml~99.9ml, step: 0.01 ml At 100.0 ml~999.9ml, step: 0.1 ml At 1000 ml~9999 ml, step: 1 ml
Time range:	00:00:01~99:59:59 (h:m:s); Step: 1s

**Micro Mode setting range:**

Parameters:	Rate, Time and Volume to Be Infused (VTBI) can be set to automatically calculate the remaining parameter.
Rate range:	0.10~100.00ml/h;
Rate steps:	0.01 ml/h
VTBI range:	0.10~1000.00ml
VTBI steps:	At 0.10 ml~99.9ml, step: 0.01 ml At 100.0 ml~999.9ml, step: 0.1 ml
Time range:	00:00:01~99:59:59 (h:m:s); Step: 1s

**Ramp mode setting range:**

Parameters: Total Time, Up Time, Down Time and Volume to be infused (VTBI) can be set to automatically calculate the Steady Rate. The rate range below refers to the values in which the automatically calculated steady rate may not exceed.

Rate range: 0.10~2000.00ml/h;

Rate steps: At 0.10 ml/h~99.9ml/h, step: 0.01 ml/h

At 100.0 ml/h~999.9ml/h, step: 0.1 ml/h

At 1000 ml/h~2000 ml/h, step: 1 ml/h

VTBI range: 0.10~9999ml

VTBI steps: At 0.10 ml~99.9ml, step: 0.01 ml

At 100.0 ml~999.9ml, step: 0.1 ml

At 1000 ml~9999 ml, step: 1 ml

Time (Total time, Up time, Down

time) range: 00:00:01~99:59:59(h: m: s); Step: 1s  
If the up time and down time are not input, it defaults to 0s.

**Intermittent Mode setting range:**

Parameters: Intmt.Vol (Intermittent Volume), Intmt. Rate (Intermittent Rate), Keep Rate, Interval and Volume to Be Infused (VTBI) can be set.

Intmt. Rate range: 0.10ml/h~2000.00ml/h

Rate steps: At 0.10 ml/h~99.9ml/h, step: 0.01 ml/h

At 100.0 ml/h~999.9ml/h, step: 0.1 ml/h

At 1000 ml/h~2000 ml/h, step: 1 ml/h

VTBI range: 0.10~9999ml

VTBI steps: At 0.10 ml~99.9ml, step: 0.01 ml

At 100.0 ml~999.9ml, step: 0.1 ml

At 1000 ml~9999 ml, step: 1 ml

Keep Rate range: 0.10~5.00ml/h; Step: 0.01ml/h

Time range: 00:00:01~99:59:59 (h:m:s); Step: 1s

**Dose Time Mode setting range:**

Parameters: Drug Volume, Drug Weight, Dose Rate, Time and Volume to be infused (VTBI) can be set to automatically calculate the Infusion rate. The rates below refer to the values in which the automatically calculated Infusion rate may not exceed.

Dose Rate range: When the unit of VTBI is other than ml, the rate is dose rate, and the range is 0.10~999 (editable). when the unit of VTBI is ml, the rate unit is ml/h, with a range of 0.10~2000.00ml/h.

Rate steps: At 0.10 ml/h~99.9ml/h, step: 0.01 ml/h

At 100.0 ml/h~999.9ml/h, step: 0.1 ml/h

At 1000 ml/h~2000 ml/h, step: 1 ml/h

VTBI range: 0.10~9999ml

VTBI steps: At 0.10 ml~99.9ml, step: 0.01 ml

At 100.0 ml~999.9ml, step: 0.1 ml

At 1000 ml~9999 ml, step: 1 ml

Time range: 00:00:01~99:59:59(h: m: s); Step: 1s

**Bolus and Purging ranges:**

Rate range: 0.10~2000.00ml/h; Step: 0.01ml/h

Volume range: 0.10~99.99ml

**Infusion Accuracy:**

Infusion accuracy:  $\leq \pm 2\%$  for rates  $> 1.00\text{ml/h}$

(For specified and calibrated syringes, accuracy of infusion:  $\leq \pm 1.8\%$ )

Bolus accuracy:  $\leq \pm 5\%$  or  $\leq \pm 0.05\text{ml}$ , whichever value is larger

**Drug Library and Information Management:**

**Drug Library:** Up to 3000 drug types can be registered.

**Patient Information Management:** A patient information management system is inbuilt that allows for patient information input and recording.

Patient Type: Adult, Neonatal, Paediatric selectable

Other information: Name, Patient ID, Bed No., DOB, Height, Weight, Age, Sex, Blood type, BSA, BMI, Doctor in Charge, Prescription

**Alarms:**

Compliance: IEC60601-1-8

Alarm levels: P1 to P4

Alarm volume: Max: 65db (At 1m distance)

Min: 45db (At 1m distance)

**Occlusion Alarm:**

Range: 150~900mmHg (P1-P4)

Accuracy:  $\pm 20\%$  or  $\pm 125\text{mmHg}$ , whichever is greater

Unit settings: mmHg, kPa, bar, psi

EBIS: 3 Levels Selectable Empty Bottle Inspection Sensitivity (EBIS)

**Features**

Syringe selection: Brands selectable upon Infusion Set installation: JieRui, BD  
Supports up to 3 user-defined custom Syringe Brands.

Syringe Parameters: 5ml, 10ml, 20ml, 30ml, 50/60ml

Accumulated

Volume Range: 0.10~9999ml

Drug Weight (Dose) 0.100~99999 (g/mg/ug/ng/IU/U/KU)

Range: Steps of 0.001

Dose Rate Range: 0.1~9999 (ng/kg/h, ug/kg/h, mg/kg/h, IU/kg/h, IU/kg/min, ug/kg/min, mg/kg/min, ng/kg/min, g/kg/h, g/kg/min, KU/kg/min, KU/kg/h, U/kg/min, U/kg/h)

Ingress Protection: IP44

KVO Rate range (Keep Vein Open):	0.10~5.00ml/h; (KVO will last for 30 minutes, then an alarm will sound) Step: 0.01ml/h	Nurse Call:	Supported
Event Logs:	Up to 2000 records FIFO. Records any system changes, alarms and actions taken. Supports record log export.	No operation time settings:	1~5min (step: 1min)
Night mode:	Automatic adjustment of volume and screen brightness	Continuous Titration:	Allows editing of infusion rate for ongoing infusions without the need to stop any currently ongoing infusions.
Purging:	Manual and automatic selectable	Others:	Automatic Pressure Release (Anti-bolus), Syringe Calibration function, WIFI
Software Update:	Via RS232 port		

**\*Notice: Specifications subject to change without prior notice. All rights reserved by COMEN**